Instructions for Use

MRX Specialty Abnormal Control, K5022

NORDIC
BIOMARKER

Contents

English [EN]	3
Dansk [DA]	6
Deutsch [DE]	9
Español [ES]	12
Magyar [HU]	15
Italiano [IT]	18
Norsk [NO]	21

Instructions for Use [EN]

MRX Specialty Abnormal Control

REF K5022

For *In vitro* Diagnostic Use.

1 Intended use

For quality control of assays intended to measure:

- PT INR
- APTT
- Fibrinogen
- Antithrombin
- D-dimer
- Free Protein S
- Protein C
- vWF Antigen

in citrated human plasma samples. Expected values are determined with MRX assays provided by Nordic Biomarker. Intended to be used by professional laboratory personnel using coagulation analysers or analysers with turbidimetric detection.

2 Components

MRX Specialty Abnormal Control consists of 10 x 1 mL:

- Lyophilised citrated human plasma.
- Buffered saline.
- Additives of protein components from bovine plasma.
- Stabiliser (BSA).
- Preservative (NaN₃).

3 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

The control contains material of human origin. Each donor has been tested by approved methods and found negative for the presence of HBsAg and anti-HIV I & II and anti-HCV. However, as no method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious.

The control contains Bovine Serum Albumin (< 5.0%) for stability, and bovine plasma. The animals were approved by veterinarians by ante- and post-mortem inspections.

However, as no method can offer complete assurance, this material should be handled as potentially infectious.

The control contains sodium azide (< 0.1%) to prevent microbial growth; use proper disposal procedures.

4 Expected values

Expected values for each new lot of MRX Specialty Abnormal Control are determined with MRX reagents on instruments as specified below:

Parameter	Reagent	Instrument
PT Owren INR	MRX PT Owren	ACL TOP
PT Quick INR	MRX PT Quick	ACL TOP
APTT	MRX APTT	ACL TOP
Fibrinogen	MRX Fib Clauss	ACL TOP
Antithrombin	MRX Antithrombin	ACL TOP
D-dimer DDU	MRX Red D-dimer	Sysmex CS-series
D-dimer FEU	MRX Red D-dimer	Sysmex CS-series
Free Protein S	MRX Blue Free Protein S	Sysmex CS-series
Protein C	MRX Protein C	Sysmex CS-series
vWF Antigen	MRX Green vWF Antigen	Sysmex CS-series

Refer to the Certificate of Analysis for the lot-specific expected values for each parameter. Expected values should be used as guidance only in the internal quality control procedure employed by the laboratory. Each laboratory is recommended to establish its own acceptable range.

5 Preparation

- Before opening, carefully tap the vial against a surface to collect the lyophilised material at the bottom.
- Add 1.00 mL deionised water (e.g. MRX Laboratory Water, K5036). The water temperature should be 15 - 25 °C.
- Reseal the vial and let it stand for approximately 15 minutes at 15 - 25 °C.
- Gently mix by swirling or rotating until the content is completely reconstituted.

6 Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 24 hours at 2 - 25 °C in the closed original vial, provided no contamination occurs.

7 Material required but not provided

- Coagulation analyser/analyser with turbidimetric detection and pipettes.
- Deionised water for reconstitution.
- Assay(s) intended to measure the specific parameter(s) in citrated human plasma samples.
- Control material with parameter values in different levels than MRX Specialty Abnormal Control.

Recommended materials are presented below:

Solution for reconstitution	REF
MRX Laboratory Water	K5036

Reagents	REF
MRX PT Owren	K5026 K5027 K5028 K5064
MRX PT Quick	K5024
MRX APTT	K5029 K5030
MRX Antithrombin	K5033
MRX Fib Clauss	K5031
MRX Red D-dimer	K5034 K5073
MRX Green D-dimer	K5011
MRX Blue D-dimer	K5035
MRX Red Free Protein S	K5009
MRX Blue Free Protein S	K5010
MRX Protein C	K5020
MRX Green vWF Antigen	K5052

Control material	REF
MRX Specialty Normal Control	K5021

8 Quality control

To maintain consistent assay results, it is recommended that controls in different levels are analysed together at regular intervals. Each laboratory is recommended to establish its own acceptable range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. A new acceptable range must be determined for each individual lot of controls. Recalibration is suggested at least when

controls are not within the acceptable range and/or each time a new lot of reagent is used.

9 Limitations

If other assays than those provided by Nordic Biomarker are employed, it must be ensured that the usage does not conflict with the intended use of the control. Refer to the corresponding reagent's instruction for use.

Non-calibrated assays (APTT assays) are particularly sensitive to variations between laboratories. Expected values should be used as guidance only.

10 Summary of Safety and Performance

The Summary of Safety and Performance (SSP) is available in the European database on medical devices (Eudamed), where it is linked to this Basic UDI-DI, 7350060320206Q. The Eudamed public website is found on <https://ec.europa.eu/tools/eudamed>. In case Eudamed is not fully functional, the SSP is available from Nordic Biomarker on request.

11 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

12 Additional information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

13 Definition of symbols

Manufacturer

CE mark

In vitro diagnostic medical device

Catalogue number

Batch code

Consult electronic instructions for use

nordicbiomarker.com/IFU

Use-by date

Temperature limit

Biological risks

Contains biological material of animal origin

Contains human blood or plasma derivatives

14 Revision history

Version	Changes to previous version
7.0	The document has been revised in connection with CE marking under Regulation (EU) 2017/746.

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