Instructions for Use

MRX PT Owren, K5064

Contents

English [EN]	3
Dansk [DA]	8
Eesti [ET]	13
Suomi [FI]	18
Lietuvių [LT]	23
Norsk [NO]	28
Svenska [SV]	33

Instructions for Use [EN]

MRX PT Owren

REF K5064

For *In vitro* Diagnostic Use.

1 Intended use

Clotting test for the quantitative determination of prothrombin time (PT) in citrated human plasma according to Owren's method. Can be used to monitor warfarin therapy and to evaluate the function of the extrinsic pathway of the coagulation cascade. Intended to be used by professional laboratory personnel using coagulation analysers.

2 Valid combinations

This device is intended for use in combination with a coagulation analyser. Valid combinations are:

Reagent	Instrument	Instrument settings
K5064	Sysmex CS-series	Refer to your local distributor for
	Sysmex CN-series	instrument-specific application sheets

Note that if the device is used in combination with other instruments or other instrument-specific application sheets not provided by Nordic Biomarker, the combination must be validated in accordance with Regulation (EU) 2017/746.

3 Background and principle of method

Prothrombin time determination using MRX PT Owren is dependent on the activity of vitamin K-dependent coagulation factors, FII, FVII and FX. Owren's PT method is based on stimulation of the extrinsic coagulation pathway via activation of FVII by thromboplastin in the presence of calcium. The activated complex (FVIIa-tissue factor) activates FX. Subsequently, FXa, in complex with FV, activates prothrombin to thrombin leading to the conversion of fibrinogen to fibrin, detected as a clot. In patients treated with vitamin K-antagonists, the activity of vitamin K-dependent coagulation factors will be lower and result in a prolonged PT.¹

When analysing PT using MRX PT Owren, sample diluted in MRX PT Buffer and reagent are mixed, and the clotting time is measured. The final dilution of sample is 1:21. The

reagent is enriched with FV and fibrinogen from bovine plasma (deficient of FII, FVII and FX), making it insensitive to variations in patient FV and fibrinogen levels. Also, due to the 1:21 sample dilution, the Owren PT method is relatively insensitive to interferences.

4 Components

MRX PT Owren consists of:

10 × 5 mL lyophilised rabbit brain thromboplastin, bovine plasma, stabilisers, and preservatives.

5 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

MRX PT Owren contains Bovine Serum Albumin and bovine plasma. The animals were approved by veterinarians by ante- and post-mortem inspections.

However, as no method can offer complete assurance, this material should be handled as potentially infectious.

MRX PT Owren contains sodium azide (less than 0.1%) to prevent microbial growth; use proper disposal procedures.

6 Preparation

MRX PT Owren should be reconstituted and allowed to mature for 24 hours prior to use (see section 7 Storage and stability). If the reagent is kept in the instrument for several days, it is important to mix the reagent at least daily before use. Continuous stirring is not necessary.

6.1 Reconstitution according to 2X method

- Allow the reagent and deionised water to reach room temperature (15 - 25 °C).
- Add 5 mL deionised water (e.g., MRX Laboratory Water, K5036) to the reagent vial.
- Keep the reagent for 10 - 15 minutes at room temperature (15 - 25 °C). It is important to mix gently by swirling at the beginning, in the middle and at the

end of this period. The reagent will dissolve into a slightly opaque colourless liquid.

- Allow the reconstituted reagent to mature 24 h at 2 - 8 °C before use.

7 Storage and stability

Store at 2 - 8 °C. Stable for 6 days (2X method) following reconstitution and 24 h incubation (Figure 1). For information regarding on-board stability, refer to the instrument-specific application sheet.

Stability data is valid provided that reconstituted reagent is stored closed in the original vial and that no contamination occurs.


Figure 1. Stability and storage temperatures for MRX PT Owren reconstituted with the 2X method.

8 Material required but not provided

- Coagulation analyser (refer to section 2 for valid instruments), and pipettes.
- Calibrator.
- Control materials with different INR for quality control.
- MRX PT Buffer for sample dilution.
- Calcium chloride (25 mM).
- Deionised water for reconstitution.

Solution for sample dilution	REF
MRX PT Buffer	K5065 K5038

Recommended materials are presented below:

Calibrators	REF
Refer to section 10.2	-

Control material	REF
MRX Routine Normal Control and MRX Routine Abnormal Control or	K5039 K5040
MRX Specialty Normal Control and MRX Specialty Abnormal Control	K5021 K5022

Other solutions	REF
Calcium chloride (25 mM) e.g., MRX Calcium Chloride	K5049 K5050 K5051
Deionised water for reconstitution e.g., MRX Laboratory Water	K5036

9 Specimen collection and preparation

Venous blood is collected in 3.2% sodium citrate at a ratio of 9 parts blood to 1 part anticoagulant (1:10 ratio). The ratio is critical. Trauma or stasis during blood sampling should be avoided. Blood should not be collected through a heparin lock or other heparinised line. Invert immediately after sampling. The presence of any clots in a specimen is a cause for rejection. Centrifuge to produce platelet-poor plasma and use for analysis. Refer to CLSI guideline H21 for further instructions on specimen collection, handling, and storage.²

10 Procedure

For detailed information, refer to the instrument operator's manual and to the instrument-specific application sheet.

10.1 Application settings

Enter the product specific settings in the instrument according to the instrument-specific application sheet.

10.2 Calibration

Each lot of MRX PT Owren must be calibrated prior to determination of INR in patients' plasma samples. During calibration, reagent lot and instrument specific values for normal clotting time (NCT) and International Sensitivity Index (ISI) required for INR determination (see Section 11 Results) are established.

Local ISI calibration can be performed using lyophilised calibrant plasmas,^{3, 4} such as the calibration kit provided by Equalis. Refer to instructions provided by the calibrator supplier for calibration procedure.

It is also possible to determine lot and instrument specific mean normal prothrombin time (MNPT) and ISI using guidelines provided by WHO.⁵

10.3 Quality control

To maintain consistent assay results, it is recommended that controls in different levels are analysed together at regular intervals. Each laboratory is recommended to establish its own acceptable range to determine the allowable variation in the day-to-day performance of the

test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. A new acceptable range must be determined for each individual lot of controls. Refer to section 8 "Material required but not provided" for recommended control materials.

10.4 Sample analysis

For analysis of samples, refer to the instrument operator's manual.

11 Results

The results are reported in International Normalised Ratio (INR). INR is calculated by dividing the obtained clotting time by the normal clotting time or mean normal PT (see Section 10.2 Calibration) raised to the lot specific ISI value for the reagent and measurement system according to the formula below.^{6,7}

$$\text{INR} = (\text{patient's PT} / \text{control PT})^{\text{ISI}}$$

Patient's PT= The patient's PT in seconds.

Control PT = The normal PT time (NCT or MNPT) in seconds. Lot specific value for the reagent and instrument combination determined during calibration.

ISI = Lot specific value for the reagent and instrument combination determined during calibration.

Note: The ISI and NCT/MNPT values are specific for each lot of reagent and each reagent/instrument combination.

12 Expected values

Expected values for normal population are INR 0.9 - 1.2⁸ while expected values for affected population (e.g. patients treated with vitamin K-antagonists) may vary depending on indication.⁹

13 Limitations and interfering substances

Results may be affected by insufficient blood sampling with a shifted ratio of sodium citrate to patient plasma. The PT results may also be affected by heparin. MRX PT Owren is affected by direct/novel oral anticoagulants (DOACs/NOACs). Therefore, analysis of plasma samples containing DOACs is not recommended using MRX PT Owren.

MRX PT Owren is insensitive to the following substances on Sysmex CS-series instruments:

Interfering substance	Tolerance
Bilirubin	Up to 40 mg/dL
Haemoglobin	Up to 1000 mg/dL
Triglycerides	Up to 1000 mg/dL
Unfractionated heparin	Up to 50 U/dL
Low molecular weight heparin	Up to 100 U/dL

Each laboratory is recommended to determine its own sensitivity to heparins.

For interfering substances for other valid instruments in section 2, refer to the instrument-specific application sheet.

14 Analytical performance characteristics

Analytical performance characteristics for MRX PT Owren on Sysmex CS-series instruments:

Repeatability:

Sample	Mean INR	CV %
Level 1	1.0	1.0
Level 2	2.7	1.4
Level 3	3.9	1.0
Level 4	6.0	1.1

The overall reproducibility CV % for MRX PT Owren is ≤ 5.5% including lot, instrument, and laboratory variability factors.¹⁰

For analytical performance characteristics for other valid instruments in section 2, refer to the instrument-specific application sheet.

15 Summary of Safety and Performance

The Summary of Safety and Performance (SSP) is available in the European database on medical devices (Eudamed), where it is linked to this Basic UDI-DI, 7350060320036Q. The Eudamed public website is found on <https://ec.europa.eu/tools/eudamed>. In case Eudamed is not fully functional, the SSP is available from Nordic Biomarker on request.

16 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

17 Additional information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

The instrument-specific application sheet is available from your local distributor.

18 References

1. HORSTI, Juha. Comparison of Quick and Owren prothrombin time with regard to the harmonisation of the International Normalised Ratio (INR) system. 2002.
2. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.
3. LINDAHL, Tomas L., et al. INR calibration of Owren-type prothrombin time based on the relationship between PT% and INR utilizing normal plasma samples. *Thrombosis and haemostasis*, 2004, 91.06: 1223-1231.
4. HILLARP, Andreas, et al. Local INR calibration of the Owren type prothrombin assay greatly improves the intra- and interlaboratory variation. *Thrombosis and haemostasis*, 2004, 91.02: 300-307.
5. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. MEETING; WORLD HEALTH ORGANIZATION; WORLD HEALTH ORGANIZATION. EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. *Who Expert Committee on Biological Standardization: Sixty-second Report*. World Health Organization, 2013.
6. LOELINGER, E. A. ICSH/ICTH recommendation for reporting prothrombin time in oral anticoagulant control. *Thromb Haemostas*, 1985, 53: 155-156.
7. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION, et al. *WHO Expert Committee on Biological Standardization: forty-eighth report*. World Health Organization, 1999.
8. NILSSON, Caroline U.; STRANDBERG, Karin; REINSTRUP, Peter. Warfarin monitoring with viscoelastic haemostatic assays, thrombin generation, coagulation factors and correlations to Owren and Quick prothrombin time. *Scandinavian Journal of Clinical and Laboratory Investigation*, 2018, 78.5: 358-364.
9. KEELING, David, et al. Guidelines on oral anticoagulation with warfarin. *Br J Haematol*, 2011, 154.3: 311-24.
10. Meijer, Piet, et al. International normalized ratio (INR) testing in Europe: between-laboratory comparability of test results obtained by Quick and Owren reagents. *Clin Chem Lab Med*. 2018 Sep 25;56(10):1698-1703.

19 Definition of symbols

Manufacturer

Consult electronic instructions for use

nordicbiomarker.com/IFU

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Biological risks

Batch code

Contains biological material of animal origin

20 Revision history

Version	Changes to previous version
3.0	The document has been revised in connection with CE marking under Regulation (EU) 2017/746.

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18 Viited

1. HORSTI, Juha. Comparison of Quick and Owren prothrombin time with regard to the harmonisation of the International Normalised Ratio (INR) system. 2002.
2. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.
3. LINDAHL, Tomas L., et al. INR calibration of Owren-type prothrombin time based on the relationship between PT% and INR utilizing normal plasma samples. *Thrombosis and haemostasis*, 2004, 91.06: 1223-1231.
4. HILLARP, Andreas, et al. Local INR calibration of the Owren type prothrombin assay greatly improves the intra-and interlaboratory variation. *Thrombosis and haemostasis*, 2004, 91.02: 300-307.
5. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. MEETING; WORLD HEALTH ORGANIZATION; WORLD HEALTH ORGANIZATION. EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. *Who Expert Committee on Biological Standardization: Sixty-second Report*. World Health Organization, 2013.
6. LOELINGER, E. A. ICSH/ICTH recommendation for reporting prothrombin time in oral anticoagulant control. *Thromb Haemostas*, 1985, 53: 155-156.
7. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION, et al. *WHO Expert Committee on Biological Standardization: forty-eighth report*. World Health Organization, 1999.

8. NILSSON, Caroline U.; STRANDBERG, Karin; REINSTRUP, Peter. Warfarin monitoring with viscoelastic haemostatic assays, thrombin generation, coagulation factors and correlations to Owren and Quick prothrombin time. *Scandinavian Journal of Clinical and Laboratory Investigation*, 2018, 78.5: 358-364.
9. KEELING, David, et al. Guidelines on oral anticoagulation with warfarin. *Br J Haematol*, 2011, 154.3: 311-24.
10. Meijer, Piet, et al. International normalized ratio (INR) testing in Europe: between-laboratory comparability of test results obtained by Quick and Owren reagents. *Clin Chem Lab Med*. 2018 Sep 25;56(10):1698-1703.

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18 Viitteet

1. HORSTI, Juha. Comparison of Quick and Owren prothrombin time with regard to the harmonisation of the International Normalised Ratio (INR) system. 2002.
2. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.
3. LINDAHL, Tomas L., et al. INR calibration of Owren-type prothrombin time based on the relationship between PT% and INR utilizing normal plasma samples. *Thrombosis and haemostasis*, 2004, 91.06: 1223-1231.
4. HILLARP, Andreas, et al. Local INR calibration of the Owren type prothrombin assay greatly improves the intra-and interlaboratory variation. *Thrombosis and haemostasis*, 2004, 91.02: 300-307.
5. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. MEETING; WORLD HEALTH ORGANIZATION; WORLD HEALTH ORGANIZATION. EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. *Who Expert Committee on Biological Standardization: Sixty-second Report*. World Health Organization, 2013.
6. LOELINGER, E. A. ICSH/ICTH recommendation for reporting prothrombin time in oral anticoagulant control. *Thromb Haemostas*, 1985, 53: 155-156.
7. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION, et al. *WHO Expert Committee on Biological Standardization: forty-eighth report*. World Health Organization, 1999.
8. NILSSON, Caroline U.; STRANDBERG, Karin; REINSTRUP, Peter. Warfarin monitoring with viscoelastic haemostatic assays, thrombin generation,

coagulation factors and correlations to Owren and Quick prothrombin time. *Scandinavian Journal of Clinical and Laboratory Investigation*, 2018, 78.5: 358-364.

9. KEELING, David, et al. Guidelines on oral anticoagulation with warfarin. *Br J Haematol*, 2011, 154.3: 311-24.
10. Meijer, Piet, et al. International normalized ratio (INR) testing in Europe: between-laboratory comparability of test results obtained by Quick and Owren reagents. *Clin Chem Lab Med*. 2018 Sep 25;56(10):1698-1703.

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18 Bibliografija

1. HORSTI, Juha. Comparison of Quick and Owren prothrombin time with regard to the harmonisation of the International Normalised Ratio (INR) system. 2002.
2. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.
3. LINDAHL, Tomas L., et al. INR calibration of Owren-type prothrombin time based on the relationship between PT% and INR utilizing normal plasma samples. *Thrombosis and haemostasis*, 2004, 91.06: 1223-1231.
4. HILLARP, Andreas, et al. Local INR calibration of the Owren type prothrombin assay greatly improves the intra-and interlaboratory variation. *Thrombosis and haemostasis*, 2004, 91.02: 300-307.
5. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. MEETING; WORLD HEALTH ORGANIZATION; WORLD HEALTH ORGANIZATION. EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. *Who Expert Committee on Biological Standardization: Sixty-second Report*. World Health Organization, 2013.
6. LOELINGER, E. A. ICSH/ICTH recommendation for reporting prothrombin time in oral anticoagulant control. *Thromb Haemostas*, 1985, 53: 155-156.
7. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION, et al. *WHO Expert Committee on Biological Standardization: forty-eighth report*. World Health Organization, 1999.
8. NILSSON, Caroline U.; STRANDBERG, Karin; REINSTRUP, Peter. Warfarin monitoring with

viscoelastic haemostatic assays, thrombin generation, coagulation factors and correlations to Owren and Quick prothrombin time. *Scandinavian Journal of Clinical and Laboratory Investigation*, 2018, 78.5: 358-364.

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18 Referanser

1. HORSTI, Juha. Comparison of Quick and Owren prothrombin time with regard to the harmonisation of the International Normalised Ratio (INR) system. 2002.
2. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.
3. LINDAHL, Tomas L., et al. INR calibration of Owren-type prothrombin time based on the relationship between PT% and INR utilizing normal plasma samples. *Thrombosis and haemostasis*, 2004, 91.06: 1223-1231.
4. HILLARP, Andreas, et al. Local INR calibration of the Owren type prothrombin assay greatly improves the intra-and interlaboratory variation. *Thrombosis and haemostasis*, 2004, 91.02: 300-307.
5. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. MEETING; WORLD HEALTH ORGANIZATION; WORLD HEALTH ORGANIZATION. EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. *Who Expert Committee on Biological Standardization: Sixty-second Report*. World Health Organization, 2013.
6. LOELINGER, E. A. ICSH/ICTH recommendation for reporting prothrombin time in oral anticoagulant control. *Thromb Haemostas*, 1985, 53: 155-156.
7. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION, et al. *WHO Expert Committee on Biological Standardization: forty-eighth report*. World Health Organization, 1999.
8. NILSSON, Caroline U.; STRANDBERG, Karin; REINSTRUP, Peter. Warfarin monitoring with viscoelastic haemostatic assays, thrombin generation, coagulation factors and correlations to Owren and Quick prothrombin time. *Scandinavian Journal of Clinical and Laboratory Investigation*, 2018, 78.5: 358-364.
9. KEELING, David, et al. Guidelines on oral anticoagulation with warfarin. *Br J Haematol*, 2011, 154.3: 311-24.
10. Meijer, Piet, et al. International normalized ratio (INR) testing in Europe: between-laboratory comparability of test results obtained by Quick and Owren reagents. *Clin Chem Lab Med*. 2018 Sep 25;56(10):1698-1703.

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10. Meijer, Piet, et al. International normalized ratio (INR) testing in Europe: between-laboratory comparability of test results obtained by Quick and Owren reagents. *Clin Chem Lab Med*. 2018 Sep 25;56(10):1698-1703.

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18 Referenser

1. HORSTI, Juha. Comparison of Quick and Owren prothrombin time with regard to the harmonisation of the International Normalised Ratio (INR) system. 2002.
2. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.
3. LINDAHL, Tomas L., et al. INR calibration of Owren-type prothrombin time based on the relationship between PT% and INR utilizing normal plasma samples. *Thrombosis and haemostasis*, 2004, 91.06: 1223-1231.
4. HILLARP, Andreas, et al. Local INR calibration of the Owren type prothrombin assay greatly improves the intra- and interlaboratory variation. *Thrombosis and haemostasis*, 2004, 91.02: 300-307.
5. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. MEETING; WORLD HEALTH ORGANIZATION; WORLD HEALTH ORGANIZATION. EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. *Who Expert Committee on Biological Standardization: Sixty-second Report*. World Health Organization, 2013.
6. LOELINGER, E. A. ICSH/ICTH recommendation for reporting prothrombin time in oral anticoagulant control. *Thromb Haemostas*, 1985, 53: 155-156.
7. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION, et al. *WHO Expert Committee on Biological Standardization: forty-eighth report*. World Health Organization, 1999.
8. NILSSON, Caroline U.; STRANDBERG, Karin; REINSTRUP, Peter. Warfarin monitoring with viscoelastic haemostatic assays, thrombin generation, coagulation factors and correlations to Owren and Quick prothrombin time. *Scandinavian Journal of Clinical and Laboratory Investigation*, 2018, 78.5: 358-364.
9. KEELING, David, et al. Guidelines on oral anticoagulation with warfarin. *Br J Haematol*, 2011, 154.3: 311-24.

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