Instructions for Use [EN]

Control Level 2 SC107K

For *In vitro* Diagnostic Use.

1. Intended use

For quality control of the assay D-Dimer Latex LRT. Can also be used for quality control of assays intended to measure PT INR, APTT, Fibrinogen and Antithrombin in citrated human plasma samples. Intended to be used by professional laboratory personnel using coagulation analysers or analysers with turbidimetric detection.

2. Components

Control Level 2 consists of 10 x 1 mL:

- Lyophilised citrated human plasma
- Buffered saline.
- Additives of protein components from bovine plasma.
- Stabiliser (BSA).
- Preservative (NaN₃)

3. Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

The control contains material of human origin. Each donor has been tested by approved methods and found negative for the presence of HBsAg and anti-HIV I & II and anti-HCV. However, as no method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious.

The control contains bovine serum albumin (< 5.0%) for stability, and bovine plasma. The animals were approved by veterinarians by ante- and post-mortem inspections. However, as no method can offer complete assurance, this material should be handled as potentially infectious.

The control contains sodium azide (< 0.1%) to prevent microbial growth; use proper disposal procedures.

4. Expected values

Values for the D-Dimer parameter in each new lot of Control Level 2 are determined using the D-Dimer Latex LRT reagent. All other assay values are determined using the manufacturer's reagents. Refer to the Certificate of Analysis for the lot-specific values. Values should be used as guidance only in the internal quality control procedure employed by the laboratory. Each laboratory is recommended to establish its own acceptable range.

5. Preparation

- Before opening, carefully tap the vial against a surface to collect the lyophilised material at the bottom.
- Add 1.00 mL deionised water. The water temperature should be 15 - 25 °C.
- Reseal the vial and let it stand for approximately 15 minutes at 15 - 25 °C.
- Gently mix by swirling or rotating until the content is completely reconstituted.

6. Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 24 hours at 2 - 25 °C in the closed original vial, provided no contamination occurs.

7. Material required but not provided

D-Dimer assay:

- Analyser with turbidimetric detection (refer to the Instructions for Use associated with the D-Dimer Latex LRT assay for valid instruments).
- The D-Dimer assay intended to be quality controlled.
- Control material with D-Dimer values in different levels than Control Level 2.
- Deionised water for reconstitution.

Reagent	REF
D-Dimer Latex LRT	120011

Control material	REF
Control Level 1	SC106K

Other assays:

- Coagulation analyser/analyser with turbidimetric detection.
- Assay(s) intended to measure the specific parameter(s) in citrated human plasma samples.
- Control material with parameter values in different levels than Control Level 2.
- Deionised water for reconstitution.

8. Quality control

To maintain consistent assay results, it is recommended that controls in different levels are analysed together at regular intervals. Each laboratory is recommended to establish its own acceptable range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. A new acceptable range must be determined for each individual lot of controls. Recalibration is suggested at least when controls are not within the acceptable range and/or each time a new lot of reagent is used.

9. Limitations

If other assays than those provided by the manufacturer are employed, it must be ensured that the usage does not conflict with the intended use of the control. Refer to the corresponding reagent's instruction for use.

Non-calibrated assays (APTT assays) are particularly sensitive to variations between laboratories. Expected values should be used as guidance only.

10. Summary of Safety and Performance

The Summary of Safety and Performance (SSP) is available in the European database on medical devices (Eudamed), where it is linked to this Basic UDI-DI, 73500603201875. The Eudamed public website is found on <https://ec.europa.eu/tools/eudamed>. In case Eudamed is not fully functional, the SSP is available from the manufacturer on request.

11. Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to the manufacturer as well as the national competent authority in which the user is established.

12. Additional information

The e-IFU (other languages) and the safety data sheet are available at www.hyphen-biomed.com. A paper copy of these Instructions for Use is available on request. Contact your local distributor.

For customer support and instrument-specific application sheet, please contact your local distributor.

13. Definition of symbols

Manufacturer

www.hyphen-biomed.com

Consult electronic instructions for use

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Biological risks

Batch code

Contains biological material of animal origin

Unique device identifier

Contains human blood or plasma derivatives

14. Revision history

Version	Changes to previous version
4.0	Added translation to Danish, Finnish, Norwegian and Swedish in electronic instructions for use.

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