Instructions for Use [EN]

D-Dimer Calibrator REF SC105K

For *In vitro* Diagnostic Use.

1. Intended use

For calibration of the D-Dimer Latex LRT reagent in ng/mL Fibrinogen Equivalent Unit (FEU). Intended to be used by professional laboratory personnel using analysers with turbidimetric detection.

2. Background and principle of method

For background information and principle of the method, refer to the Instructions for Use associated with the D-Dimer Latex LRT assay.

3. Components

D-Dimer Calibrator consists of:
4 x 1 mL lyophilised citrated human plasma enriched with D-Dimer.

4. Metrological Traceability

The metrological traceability for D-Dimer Calibrator is established in accordance with EN ISO 17511:2021.¹ An in-house reference material for D-Dimer Calibrator is established since no international conventional reference material is available for D-Dimer. D-Dimer Calibrator fulfils the requirement for maximum allowable expanded measurement uncertainty based on biological variation.

5. Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

The calibrator contains material of human origin. Each donor has been tested by approved methods and found negative for the presence of HBsAg and anti-HIV I & II and anti-HCV. However, as no method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious.

6. Preparation

- Before opening, carefully tap the vial against a surface to collect the lyophilised material at the bottom.
- Add 1.00 mL deionised water. The water temperature should be 15 - 25 °C.
- Reseal the vial and let it stand for approximately 15 minutes at 15 - 25 °C.
- Gently mix by swirling or rotating until the content is completely reconstituted.

7. Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 24 h at 2 - 25 °C in the closed original vial, provided no contamination occurs.

8. Material required but not provided

- Analyser with turbidimetric detection (refer to the Instructions for Use associated with the D-Dimer Latex LRT assay for valid instruments) and pipettes.
- The D-Dimer reagent intended to be calibrated (D-Dimer Latex LRT).
- Control materials with different D-Dimer levels for quality control.
- Deionised water for reconstitution of calibrator and control materials.
- Phosphate buffered saline (PBS) for dilution.

Reagent	REF
D-Dimer Latex LRT	120011

Recommended materials are presented below:

Control material	REF
Control Level 1	SC106K
Control Level 2	SC107K

Solution for dilution	REF
D-Dimer PBS Diluent	AR037K

9. Procedure

For detailed information, refer to the instrument operator's manual and to the instrument-specific application sheet.

9.1 Application settings

Enter the product specific settings in the instrument according to the instrument-specific application sheet associated with the D-Dimer Latex LRT assay.

9.2 Calibration

Each lot of D-Dimer Latex LRT must be calibrated prior to determination of D-Dimer values in patients' plasma samples. Calibrate the assay according to the instrument-specific application sheet.

The lot-specific D-Dimer concentration of the D-Dimer Calibrator is found in the Certificate of Analysis.

Specification:

Analyte	Specification
D-Dimer FEU	6750 – 7750 ng/mL

Recalibration must be performed when a new lot of D-Dimer Latex LRT is introduced, and it is also suggested when controls are not within the acceptable range.

9.3 Quality control

To maintain consistent assay results, it is recommended that controls in different levels are analysed together at regular intervals. Each laboratory is recommended to establish its own acceptable range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. A new acceptable range must be determined for each individual lot of controls. Refer to section 8 "Material required but not provided" for recommended control materials.

9.4 Sample analysis

For analysis of samples, refer to the instrument operator's manual.

10. Results

The results are reported in ng/mL Fibrinogen Equivalent Unit (FEU).

Results in ng/mL can be converted to mg/L or µg/mL.

Ex. 500 ng/mL FEU = 0.5 mg/L FEU = 0.5 µg/mL FEU

11. Summary of Safety and Performance

The Summary of Safety and Performance (SSP) is available in the European database on medical devices (Eudamed), where it is linked to this Basic UDI-DI, 7350060320297A. The Eudamed public website is found on <https://ec.europa.eu/tools/eudamed>. In case Eudamed is not fully functional, the SSP is available from the manufacturer on request.

12. Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to the manufacturer as well as the national competent authority in which the user is established.

13. Additional information

The e-IFU (other languages) and the safety data sheet are available at www.hyphen-biomed.com. A paper copy of these Instructions for Use is available on request. Contact your local distributor.

For customer support and instrument-specific application sheet, please contact your local distributor.

14. References

1. EN ISO 17511:2021 In vitro diagnostic medical devices – Requirements for establishing metrological traceability of values assigned to calibrators, trueness control materials and human samples.

15. Definition of symbols

	Manufacturer		Consult electronic instructions for use
	CE mark		Use-by date
	In vitro diagnostic medical device		Temperature limit
	Catalogue number		Biological risks
	Batch code		Contains biological material of animal origin
	Unique device identifier		Contains human blood or plasma derivatives

16. Revision history

Version	Changes to previous version
4.0	Added translation to Danish, Finnish, Norwegian and Swedish in electronic instructions for use.

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