Instructions for Use

MRX Antithrombin Calibrator, K5041

Contents

English [EN]	3
Čeština [CS]	6
Dansk [DA]	9
Deutsch [DE]	12
Español [ES]	15
Soumi [FI]	18
Hrvatski [HR]	21
Magyar [HU]	24
Italiano [IT]	27
Norsk [NO]	30
Slovenčina [SK]	33
Svenska [SV]	36

Instructions for Use [EN]

MRX Antithrombin Calibrator

REF K5041

For *In vitro* Diagnostic Use.

1 Intended use

For calibration of the MRX Antithrombin reagent in IU/mL. Intended to be used by professional laboratory personnel using coagulation analysers.

2 Background and principle of method

For background information and principle of the method, refer to the Instructions for Use associated with MRX Antithrombin (K5033).

3 Components

MRX Antithrombin Calibrator consists of:
4 × 1 mL lyophilised citrated human plasma with a specified level of antithrombin.

4 Metrological traceability

Each new lot of MRX Antithrombin Calibrator is assigned with MRX Antithrombin against SSC/ISTH Secondary Coagulation Standard with traceability to WHO International Standard (NBISC code: 08/258). Refer to the Certificate of Analysis or the label for the lot-specific Antithrombin concentration.

Specification:

Analyte	Specification
Antithrombin	0.90-1.10 IU/mL

5 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

The calibrator contains material of human origin. Each donor has been tested by approved methods and found negative for the presence of HBsAg and anti-HIV I & II and anti-HCV. However, as no method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious.

6 Preparation

- Allow the vial to equilibrate at 15 - 25 °C for at least 10 minutes before reconstitution.
- Before opening, carefully tap the vial against a surface to collect the lyophilised material at the bottom.
- Add 1.00 mL deionised water (e.g., MRX Laboratory Water, K5036). The water temperature should be 15 - 25 °C.
- Reseal the vial and swirl gently. Keep reconstituted calibrator at 15 - 25 °C for at least 15 minutes before use.
- Make sure that the content is completely reconstituted and mix thoroughly by swirling before use.

7 Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 8 hours at 15 - 25 °C in the closed original vial, provided no contamination occurs.

8 Specimen collection and preparation

For specimen collection and preparation, refer to the Instructions for Use associated with MRX Antithrombin (K5033).

9 Procedure

For each instrument, refer to its operator's manual and to the instrument-specific application sheet.

10 Material required but not provided

Coagulation analyser capable of detection at 405 nm, pipettes, and the following:

Reagent	REF
MRX Antithrombin	K5033

Control material	REF
MRX Routine Normal Control	K5039
MRX Routine Abnormal Control	K5040

Solutions	REF
Deionised water for reconstitution, e.g. MRX Laboratory Water	K5036
Saline solution for dilution, e.g. MRX Sodium Chloride Diluent	K5046

17 Additional Information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

The instrument-specific application sheet is available from your local distributor.

11 Quality control

To maintain consistent assay results, it is recommended that control plasmas are assayed at regular intervals. MRX Routine Controls (K5039/K5040) are recommended for MRX Antithrombin. Each laboratory should establish a control range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. Recalibration is suggested, as a minimum, whenever control plasmas are not within the acceptable range and each time a new batch of reagent is used.

12 Results

The results are reported in IU/mL. 100% AT activity is defined as 1 IU AT function in 1 mL plasma of the WHO 3rd International Standard (NIBSC code 08/258).

13 Expected values

For expected values, refer to the Instructions for Use associated with MRX Antithrombin (K5033).

14 Limitations and interfering substances

For limitations and interfering substances, refer to the Instructions for Use associated with MRX Antithrombin (K5033).

15 Analytical performance characteristics

For analytical performance characteristics, refer to the Instructions for Use associated with MRX Antithrombin (K5033).

16 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

18 Definition of symbols

Manufacturer

Consult electronic instructions for use

nordicbiomarker.com/IFU

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Biological risks

Batch code

Contains human blood or plasma derivatives

19 Revision history

Version	Changes to previous version
4.0	Added translations to Croatian, Czech, Hungarian, Italian and Slovak.

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