Instructions for Use

MRX Routine Normal Control, K5039

Contents

| | |
|------------------|----|
| English [EN] | 3 |
| Čeština [CS] | 6 |
| Dansk [DA] | 9 |
| Deutsch [DE] | 12 |
| Español [ES] | 15 |
| Eesti [ET] | 18 |
| Soumi [FI] | 21 |
| Hrvatski [HR] | 24 |
| Magyar [HU] | 27 |
| Italiano [IT] | 30 |
| Lietuvių [LT] | 33 |
| Latviešu [LV] | 36 |
| Norsk [NO] | 39 |
| Slovenčina [SK] | 42 |
| Slovenščina [SL] | 45 |
| Svenska [SV] | 48 |

Instructions for Use [EN]

MRX Routine Normal Control

REF K5039

For *In vitro* Diagnostic Use.

1 Intended use

For quality control of assays intended to measure:

- PT INR
- Thrombin Time
- APTT
- Fibrinogen
- Antithrombin
- D-dimer

in citrated human plasma samples. Expected values are determined with MRX assays provided by Nordic Biomarker. Intended to be used by professional laboratory personnel using coagulation analysers or analysers with turbidimetric detection.

2 Components

MRX Routine Normal Control consists of 10 x 1 mL:

- Lyophilised citrated human plasma.
- Buffer.

3 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

The control contains material of human origin. Each donor has been tested by approved methods and found negative for the presence of HBsAg, anti-HIV I & II and anti-HCV. However, as no method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious.

4 Expected values

Expected values for each new lot of MRX Routine Normal Control are determined with MRX reagents on instruments as specified below:

| Parameter | Reagent | Instrument |
|---------------|-------------------|------------------|
| PT Owren INR | MRX PT Owren | ACL TOP |
| PT Quick INR | MRX PT Quick | ACL TOP |
| Thrombin Time | MRX Thrombin Time | ACL TOP |
| APTT | MRX APTT | ACL TOP |
| Fibrinogen | MRX Fib Clauss | ACL TOP |
| Antithrombin | MRX Antithrombin | ACL TOP |
| D-dimer DDU | MRX Red D-dimer | Sysmex CS-series |
| D-dimer FEU | MRX Red D-dimer | Sysmex CS-series |

Refer to the Certificate of Analysis for the lot-specific expected values for each parameter. Expected values should be used as guidance only in the internal quality control procedure employed by the laboratory. Each laboratory is recommended to establish its own acceptable range.

5 Preparation

- Before opening, carefully tap the vial against a surface to collect the lyophilised material at the bottom.
- Add 1.00 mL deionised water (e.g. MRX Laboratory Water, K5036). The water temperature should be 15 - 25 °C.
- Reseal the vial and let it stand for approximately 15 minutes at 15 - 25 °C.
- Gently mix by swirling or rotating until the content is completely reconstituted.

6 Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 24 hours at 2 - 25 °C in the closed original vial, provided no contamination occurs.

7 Material required but not provided

- Coagulation analyser/analyser with turbidimetric detection and pipettes.
- Deionised water for reconstitution.
- Assay(s) intended to measure the specific parameter(s) in citrated human plasma samples.
- Control material with parameter values in different levels than MRX Routine Normal Control.

Recommended materials are presented below:

| Solution for reconstitution | REF |
|-----------------------------|-------|
| MRX Laboratory Water | K5036 |

| Reagents | REF |
|-------------------|----------------------------------|
| MRX PT Owren | K5026 K5027 K5028 K5064 |
| MRX PT Quick | K5024 |
| MRX APTT | K5029 K5030 |
| MRX Antithrombin | K5033 |
| MRX Fib Clauss | K5031 |
| MRX Thrombin Time | K5032 |
| MRX Red D-dimer | K5034 K5073 |
| MRX Green D-dimer | K5011 |
| MRX Blue D-dimer | K5035 |

| Control material | REF |
|------------------------------|-------|
| MRX Routine Abnormal Control | K5040 |

8 Quality control

To maintain consistent assay results, it is recommended that controls in different levels are analysed together at regular intervals. Each laboratory is recommended to establish its own acceptable range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. A new acceptable range must be determined for each individual lot of controls. Recalibration is suggested at least when controls are not within the acceptable range and/or each time a new lot of reagent is used.

9 Limitations

If other assays than those provided by Nordic Biomarker are employed, it must be ensured that the usage does not conflict with the intended use of the control. Refer to the corresponding reagent's instruction for use.

Non-calibrated assays (APTT- and Thrombin Time assays) are particularly sensitive to variations between laboratories. Expected values should be used as guidance only.

10 Summary of Safety and Performance

The Summary of Safety and Performance (SSP) is available in the European database on medical devices (Eudamed), where it is linked to this Basic UDI-DI, 73500603201773. The Eudamed public website is found on <https://ec.europa.eu/tools/eudamed>. In case Eudamed is not fully functional, the SSP is available from Nordic Biomarker on request.

11 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

12 Additional information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

13 Definition of symbols

Manufacturer

Consult electronic instructions for use

nordicbiomarker.com/IFU

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Biological risks

Batch code

Contains human blood or plasma derivatives

14 Revision history

| Version | Changes to previous version |
|---------|---|
| 7.0 | Added translations to Croatian, Hungarian and Slovak. |

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