Instructions for Use

MRX Fib Clauss, K5031

Contents

English [EN]	3
Español [ES]	7

Instructions for Use [EN]

MRX Fib Clauss

REF K5031For *In vitro* Diagnostic Use.

1 Intended use

Clotting test for quantitative determination of fibrinogen in citrated human plasma according to the Clauss method. Can be used as an aid in clinical investigations when suspecting abnormal fibrinogen levels. Intended to be used by professional laboratory personnel using coagulation analysers.

2 Background and principle of method

Fibrinogen is a 340 kDa soluble plasma glycoprotein that is converted by thrombin to fibrin during blood clot formation. In unexplained bleeding cases, it can be of clinical importance to quantitate fibrinogen.¹

MRX Fib Clauss utilises the fibrinogen Clauss method where a high concentration of thrombin is added to diluted test plasma² and the clotting time is measured. In the presence of an excess of thrombin, fibrinogen is transformed into fibrin and the clot formation time is inversely proportional to the concentration of fibrinogen in the plasma sample. The test result is compared with a calibration curve prepared by clotting a series of dilutions of a reference plasma sample of known fibrinogen concentration, and a result in g/L is obtained.

3 Components

MRX Fib Clauss consists of:
10 × 5 mL lyophilised bovine thrombin (100 UNIH/mL) with Bovine Serum Albumin, buffer, stabilisers, and preservative.

4 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

MRX Fib Clauss contains sodium azide (less than 0.1%) to prevent microbial growth; use proper disposal procedures.

MRX Fib Clauss contains material of animal origin and should be handled as potentially infectious.

EUH210: Safety data sheet available on request.

5 Preparation

- Allow the vial to equilibrate at 15 - 25 °C for at least 10 minutes before reconstitution.
- Dissolve the content of each vial with 5 mL of deionised water (e.g., MRX Laboratory Water, K5036).
- Swirl gently. Make sure of the complete reconstitution of the product.
- Keep the reagent at 15 - 25 °C for 30 minutes and invert to mix before use. Do not shake.

6 Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 7 days at 2 - 25 °C in the closed original vial, provided no contamination occurs.

On board stability: 3 days at 15 °C (stored on board in open vial during daytime, 8 am - 5 pm, closed vial at 2 - 8 °C during night-time).

7 Specimen collection and preparation

Venous blood is collected in 3.2% sodium citrate at a ratio of 9 parts blood to 1 part anticoagulant (1:10 ratio). The ratio is critical. Trauma or stasis during blood sampling should be avoided. Blood should not be collected through a heparin lock or other heparinised line. Inverse immediately after sampling. The presence of any clots in a specimen is a cause for rejection. Centrifuge to produce platelet-poor plasma and use for analysis. Refer to CLSI guideline H21 for further instructions on specimen collection, handling and storage.³

8 Procedure

8.1 Calibration

Each lot of MRX Fib Clauss must be calibrated prior to determination of fibrinogen concentration in patients' plasma samples.

It is recommended to use MRX Fib Clauss Calibrator (K5042) for calibration. For calibration on automated instruments, see instrument specific application sheets available from your local distributor. For calibration on semi-automated and manual instruments, see below.

8.2 Automated instruments

For each instrument, refer to its operator's manual and to the instrument-specific application sheet.

8.3 Semi-automated and manual instruments

Make sure that the final volume in the cuvette is within instrument specification. All test tubes and pipette tips should be in plastic or in siliconised glass.

Calibration

- Using MRX Sodium Chloride Diluent (K5046), prepare dilutions of MRX Fib Clauss Calibrator: 1:5, 1:10, 1:20, 1:30 and 1:40. Perform determinations on each dilution of the MRX Fib Clauss Calibrator as follows:
- Incubate 200 µL of diluted calibrator for 2 minutes at 37 °C.
- Add 100 µL of MRX Fib Clauss reagent and start recording the time.
- Obtain the clotting times for each of the dilutions of the MRX Fib Clauss Calibrator.
- Plot the average clotting times obtained versus the respective fibrinogen concentration.

Analysis

- Dilute the test plasma 1:10 in MRX Sodium Chloride Diluent (K5046).
- Incubate 200 µL of diluted plasma sample for 2 minutes at 37 °C.
- Add 100 µL of MRX Fib Clauss reagent to the diluted plasma sample and start recording the time.
- Record the clotting time in seconds.

9 Material required but not provided

Coagulation analyser, cuvettes, pipettes, and the following:

Calibrator	REF
MRX Fib Clauss Calibrator	K5042
Control material	REF
MRX Routine Normal Control	K5039
MRX Routine Abnormal Control	K5040
Solutions	REF
Saline solution for dilution, e.g. MRX Sodium Chloride Diluent	K5046
Deionised water for reconstitution e.g. MRX Laboratory Water	K5036

10 Quality control

To maintain consistent assay results, it is recommended that control plasmas are assayed at regular intervals. MRX Routine Controls (K5039/K5040) are recommended for MRX Fib Clauss. Each laboratory should establish a control range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice.⁴

11 Results

The results are reported in g/L.

12 Expected values

Generally, the expected values are 1.5 to 3.5 g/L.⁴

Results from fibrinogen concentrations determination in plasma samples from 92 healthy blood donors using MRX Fib Clauss are presented below. The analysis was performed using an ACL Top 700 CTS instrument.

N	Mean Fibrinogen	Mean ± 2 SD
92 samples	2.53 g/L	1.49 - 3.57 g/L

Due to inter-laboratory variation, each laboratory should set its own reference interval.

13 Limitations and interfering substances

Degradation products of fibrin(ogen) in the plasma sample may affect the results. Haemolysed plasma could not be assayed.

MRX Fib Clauss is insensitive to the following substances on ACL Top 700 CTS instrument:

Interfering substance	Tolerance
Bilirubin	Up to 5 mg/dL
Triglycerides	Up to 375 mg/dL
Dabigatran	Up to 0.0475 mg/dL
Unfractionated heparin	Up to 330 U/dL

14 Analytical performance characteristics

The following performance data was obtained with an ACL Top 700 CTS instrument. Performance will depend on the instrument used.

MRX Fib Clauss has a measuring range of 0.6 - 8.0 g/L. When compared to another Fib Clauss reagent, MRX Fib Clauss correlates as follows:

y (MRX Fib Clauss on ACL Top) = $1.2 \times$ (HemosIL QFA on ACL Top) - 0.4; $r^2 = 0.98$.

Precision:

Sample	Mean Fibrinogen	Repeatability CV
Level 1	1.23 g/L	3.2%
Level 2	2.45 g/L	4.0%

15 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

16 Additional information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

The instrument-specific application sheet is available from your local distributor.

17 References

1. DAVIS, E. W.; HOUGIE, C.; LUNDBLAD, R. L. Mechanisms of blood coagulation. *recent advances in blood coagulation*. J & A Churchill Ltd, London, 1969.
2. CLAUSS, Arnold. Gennungsphysiologische schnellmethode sur Bestimmung des Fibrinogens. *Acta haematol*, 1957, 17: 237-246.
3. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.

4. Clinical and Laboratory Standards Institute. Procedure for the determination of fibrinogen in plasma; Approved Guideline, Second Edition, CLSI Document H30-A2; Vol. 21 No.18

18 Definition of symbols

Manufacturer

Consult electronic instructions for use

nordicbiomarker.com/IFU

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Biological risks

Batch code

Contains biological material of animal origin

19 Revision history

Version	Changes to previous version
4.0	Section 4: Information on material of animal origin updated. Minor editorial revisions made.

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17 Referencias

1. DAVIS, E. W.; HOUGIE, C.; LUNDBLAD, R. L. Mechanisms of blood coagulation. *recent advances in blood coagulation*. J & A Churchill Ltd, London, 1969.
2. CLAUSS, Arnold. Gennungsphysiologishe schnellmethode sur Bestimmung des Fibrinogens. *Acta haematol*, 1957, 17: 237-246.
3. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.
4. Clinical and Laboratory Standards Institute. Procedure for the determination of fibrinogen in plasma; Approved Guideline, Second Edition, CLSI Document H30-A2; Vol. 21 No.18

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