Instructions for Use

MRX Protein C, K5020

Contents

English [EN] 3

Español [ES] 7

Italiano [IT] 11

Instructions for Use [EN]

MRX Protein C

REF K5020

For *In vitro* Diagnostic Use.

1 Intended use

Chromogenic test for quantitative determination of Protein C activity in citrated human plasma as an aid in diagnosis of Protein C deficiency. Intended to be used by professional laboratory personnel using coagulation analysers measuring at 405 nm.

2 Background and principle of method

Protein C is a vitamin K dependent glycoprotein with anticoagulant properties.^{1,2} When activated by thrombin, activated Protein C (APC) regulates the coagulation process by inactivating factors Va and VIIIa. This disrupts the tenase- and prothrombinase complexes, which inhibits thrombin activation.^{3,4,5}

Protein C deficiencies, hereditary and/or acquired, increases the risk for developing thromboembolic disease.⁶ There are two different types of Protein C deficiency; type I refers to reduced biosynthesis whereas type II is characterised by a decrease in functional activity.⁷

MRX Protein C measures Protein C activity, thus covering both types of Protein C deficiency. The principle involves activation of all functional Protein C present in a plasma sample by a serine protease extracted from the venom of *Agkistrodon contortrix*. After incubation, a chromogenic substrate which is readily hydrolysed by APC is added. The pNA released is proportional to the Protein C level in the sample, which can be measured at 405 nm.

3 Components

MRX Protein C consists of:

- Activator: 3 × 2.5 mL buffer containing stabilisers and preservatives.
- Substrate: 3 × 2.5 mL chromogenic substrate containing stabilisers and preservatives.

4 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

The Activator contains material of animal origin and should be handled as potentially infectious.

The Activator and Substrate contain sodium azide (less than 0.1%) and 2-methylisothiazol-3(2H)-one (less than 0.0015%) to prevent microbial growth; use proper disposal procedures.

EUH208: Contains 2-methylisothiazol-3(2H)-one. May produce an allergic reaction.

EUH210: Safety data sheet available on request.

5 Preparation

- Activator: Ready to use. Swirl the vial gently a few times before use. Do not shake.
- Substrate: Ready to use. Swirl the vial gently a few times before use. Do not shake.

6 Storage and stability

- Activator: Store at 2 - 8 °C. Do not freeze. After opening, stable for 8 weeks at 2 - 8 °C in the closed original vial, provided no contamination occurs.
- Substrate: Store at 2 - 8 °C. Do not freeze. After opening, stable for 8 weeks at 2 - 8 °C in the closed original vial, provided no contamination occurs. Discard if the substrate solution appears yellow.

Note: Protein C substrate is light-sensitive and should be stored in the dark.

7 Specimen collection and preparation

Venous blood is collected in 3.2% sodium citrate at a ratio of 9 parts blood to 1 part anticoagulant (1:10 ratio). The ratio is critical. Trauma or stasis during blood sampling should be avoided. Inverse immediately after sampling. The presence of any clots in a specimen is a cause for rejection. Centrifuge to produce platelet-poor plasma and use for analysis. Refer to CLSI guideline H21 for further instructions on specimen collection, handling and storage.⁸

8 Procedure

For each instrument, refer to its operator's manual and to the instrument-specific application sheet.

9 Material required but not provided

Coagulation analyser capable of detection 405 nm, pipettes, and the following:

Calibrator	REF
MRX Protein C Calibrator	K5023

Control material	REF
MRX Specialty Normal Control	K5021
MRX Specialty Abnormal Control	K5022

Solutions	REF
Phosphate buffered saline (PBS) for dilution, e.g. MRX PBS Diluent	K5047
Deionised water for reconstitution e.g. MRX Laboratory Water	K5036

10 Quality control

To maintain consistent assay results, it is recommended that control plasmas are assayed at regular intervals. MRX Specialty controls (K5021/K5022) are recommended for MRX Protein C. Each laboratory should establish a control range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. Recalibration is suggested, as a minimum, whenever control plasmas are not within the acceptable range and each time a new batch of reagent is used.

11 Results

The results are reported in % Protein C, where 100% corresponds to 1 IU/mL Protein C of the WHO 2nd International Standard (NIBSC code 02/342).

12 Expected values

Normal Protein C levels are typically in the 70 - 150% range, while the levels for males are somewhat higher than for females.⁹

Results from Protein C concentrations determination in plasma samples from 120 healthy blood donors using MRX Protein C are presented below. The analysis was performed using a Sysmex CS-2100i instrument.

No of samples	2.5 th – 97.5 th percentile
120	72 - 148% PC

Due to inter-laboratory variation, each laboratory should set its own reference interval.

13 Limitations and interfering substances

The results should be used together with other clinical and diagnostic information in forming a diagnosis and for patient management.

Turbid or opalescent plasma may cause erratic results and should be interpreted with caution: dilute the sample and re-assay. MRX Protein C is insensitive to the following substances on Sysmex CS-instrument series:

Interfering substance	Tolerance
Conjugated Bilirubin	Up to 25 mg/dL
Unconjugated Bilirubin	Up to 25 mg/dL
Haemoglobin	Up to 240 mg/dL
Triglycerides	Up to 375 mg/dL
Unfractionated heparin	Up to 200 U/dL

Warfarin and aprotinin causes inhibition of activated Protein C, thus, it is expected to observe low Protein C activity in warfarin/aprotinin treated patients.

14 Analytical performance characteristics

The following performance data was obtained mainly with a Sysmex CS-2500 instrument. Correlation was tested on Sysmex CS-2100i. Reproducibility was tested on both Sysmex CS-2100i and CS-2500. Performance will depend on the instrument used.

MRX Protein C has a measuring range of 8 - 186% Protein C.

Method comparison performed using Passing-Bablok:

Predicate Device (x)	System (y)	Regression	R
Hemos IL Protein C	MRX Protein C	$y = 1.04x - 1.96$	0.99

Precision:

Sample	Repeatability		Reproducibility	
	% PC	CV %	% PC	CV %
Level 1	18	2.8	19	3.5
Level 2	38	1.5	40	3.3
Level 3	103	1.4	106	2.5

15 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

16 Additional information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.
The instrument-specific application sheet is available from your local distributor.

17 References

1. DAHLBÄCK, Björn; VILLOUTREIX, Bruno O. Regulation of blood coagulation by the protein C anticoagulant pathway: novel insights into structure-function relationships and molecular recognition. *Arteriosclerosis, thrombosis, and vascular biology*, 2005, 25.7: 1311-1320.

2. WALKER, Frederick J.; FAY, Philip J. Regulation of blood coagulation by the protein C system. *The FASEB journal*, 1992, 6.8: 2561-2567.

3. ESMON, Charles T. The protein C pathway. *Chest*, 2003, 124.3: 26S-32S.

4. DAHLBÄCK, Björn; VILLOUTREIX, Bruno O. The anticoagulant protein C pathway. *FEBS letters*, 2005, 579.15: 3310-3316.

5. STAVENUITER, Fabian; BOUWENS, Eveline AM; MOSNIER, Laurent O. Down-regulation of the clotting cascade by the protein C pathway. In: *Hematology education/... Congress of the European Hematology Association. European Hematology Association. Congress. Education Program*. NIH Public Access, 2013. p. 365.

6. GOLDENBERG, N. A.; MANCO-JOHNSON, Marilyn J. Protein C deficiency. *Haemophilia*, 2008, 14.6: 1214-1221.

7. REITSMA, P. H., et al. Protein C deficiency: a database of mutations, 1995 update. *Thrombosis and haemostasis*, 1995, 73.05: 876-889.

8. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.

9. PAGANA, Kathleen Deska; PAGANA, Timothy J. *Mosby's Manual of Diagnostic and Laboratory Tests-E-Book*. Elsevier Health Sciences, 2017.

18 Definition of symbols

	Manufacturer		Consult electronic instructions for use
		nordicbiomarker.com/IFU	
	CE mark		Use-by date
	In vitro diagnostic medical device		Temperature limit
	Catalogue number		Biological risks
	Batch code		Contains biological material of animal origin

19 Revision history

Version	Changes to previous version
8.0	Section 4: Information on material of animal origin updated. Minor editorial revisions made.

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17 Referencias

1. DAHLBÄCK, Björn; VILLOUTREIX, Bruno O. Regulation of blood coagulation by the protein C anticoagulant pathway: novel insights into structure-function relationships and molecular recognition. *Arteriosclerosis, thrombosis, and vascular biology*, 2005, 25.7: 1311-1320.
2. WALKER, Frederick J.; FAY, Philip J. Regulation of blood coagulation by the protein C system. *The FASEB journal*, 1992, 6.8: 2561-2567.
3. ESMON, Charles T. The protein C pathway. *Chest*, 2003, 124.3: 26S-32S.
4. DAHLBÄCK, Björn; VILLOUTREIX, Bruno O. The anticoagulant protein C pathway. *FEBS letters*, 2005, 579.15: 3310-3316.

5. STAVENUITER, Fabian; BOUWENS, Eveline AM; MOSNIER, Laurent O. Down-regulation of the clotting cascade by the protein C pathway. In: *Hematology education/... Congress of the European Hematology Association. European Hematology Association. Congress. Education Program*. NIH Public Access, 2013. p. 365.
6. GOLDENBERG, N. A.; MANCO-JOHNSON, Marilyn J. Protein C deficiency. *Haemophilia*, 2008, 14.6: 1214-1221.
7. REITSMA, P. H., et al. Protein C deficiency: a database of mutations, 1995 update. *Thrombosis and haemostasis*, 1995, 73.05: 876-889.
8. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.
9. PAGANA, Kathleen Deska; PAGANA, Timothy J. *Mosby's Manual of Diagnostic and Laboratory Tests-E-Book*. Elsevier Health Sciences, 2017.

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18 Bibliografia

1. DAHLBÄCK, Björn; VILLOUTREIX, Bruno O. Regulation of blood coagulation by the protein C anticoagulant pathway: novel insights into structure-function relationships and molecular recognition. *Arteriosclerosis, thrombosis, and vascular biology*, 2005, 25.7: 1311-1320.
2. WALKER, Frederick J.; FAY, Philip J. Regulation of blood coagulation by the protein C system. *The FASEB journal*, 1992, 6.8: 2561-2567.
3. ESMON, Charles T. The protein C pathway. *Chest*, 2003, 124.3: 26S-32S.
4. DAHLBÄCK, Björn; VILLOUTREIX, Bruno O. The anticoagulant protein C pathway. *FEBS letters*, 2005, 579.15: 3310-3316.
5. STAVENUITER, Fabian; BOUWENS, Eveline AM; MOSNIER, Laurent O. Down-regulation of the clotting cascade by the protein C pathway. In: *Hematology education/... Congress of the European Hematology Association. European Hematology Association. Congress. Education Program*. NIH Public Access, 2013. p. 365.
6. GOLDENBERG, N. A.; MANCO-JOHNSON, Marilyn J. Protein C deficiency. *Haemophilia*, 2008, 14.6: 1214-1221.
7. REITSMA, P. H., et al. Protein C deficiency: a database of mutations, 1995 update. *Thrombosis and haemostasis*, 1995, 73.05: 876-889.
8. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.
9. PAGANA, Kathleen Deska; PAGANA, Timothy J. *Mosby's Manual of Diagnostic and Laboratory Tests-E-Book*. Elsevier Health Sciences, 2017.

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