Instructions for Use

MRX Thrombin Time, K5032

Contents

English [EN]	3
Español [ES]	6

Instructions for Use [EN]

MRX Thrombin Time

REF K5032For *In vitro* Diagnostic Use.

1 Intended use

Clotting test for determination of the thrombin time in citrated human plasma. Intended to be used by professional laboratory personnel using coagulation analysers.

2 Background and principle of method

Thrombin time is a clotting test that measures the ability of exogenous thrombin to convert fibrinogen to fibrin, followed by fibrin polymerisation and clot detection. The time taken from addition of thrombin to citrated plasma, to fibrin clot formation, is recorded as the thrombin time and presented in seconds.

3 Components

MRX Thrombin Time consists of:

6 × 2 mL lyophilised bovine thrombin (approximately 4 UNIH/mL) Bovine Serum Albumin, buffer, stabilisers, and preservative.

4 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

MRX Thrombin Time contains sodium azide (less than 0.1%) to prevent microbial growth; use proper disposal procedures.

MRX Thrombin Time contains material of animal origin and should be handled as potentially infectious.

5 Preparation

- Allow the vial to equilibrate at 15 - 25 °C for at least 10 minutes before reconstitution.
- Reconstitute the content of each vial with 2 mL of deionised water (e.g. MRX Laboratory Water, K5036).
- Keep the reagent at 15 - 25 °C for 30 minutes.
- Visually inspect the solution, confirm that all material is dissolved. Invert several times to mix before use. Do not shake.

6 Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 7 days at 2 - 8 °C and 2 days at 15 - 25 °C in the closed original vial, provided no contamination occurs.

7 Specimen collection and preparation

Venous blood is collected in 3.2% sodium citrate at a ratio of 9 parts blood to 1 part anticoagulant (1:10 ratio). The ratio is critical. Trauma or stasis during blood sampling should be avoided. Blood should not be collected through a heparin lock or other heparinised line. Inverse immediately after sampling. The presence of any clots in a specimen is a cause for rejection. Centrifuge to produce platelet-poor plasma and use for analysis. Refer to CLSI guideline H21 for further instructions on specimen collection, handling and storage.¹

8 Procedure

8.1 Automated instruments

For each instrument, refer to its operator's manual and to the instrument-specific application sheet.

8.2 Semi-automated and manual instruments

Make sure that the final volume in the cuvette is within instrument specification. All test tubes and pipette tips should be in plastic or in siliconised glass.

- Add 100 µL plasma sample to an instrument cuvette.
- Incubate for 2 minutes at 37 °C.
- Add 50 µL MRX Thrombin Time reagent to the plasma sample and start recording the time.
- Record the clotting time in seconds.

9 Material required but not provided

Coagulation analyser, cuvettes, pipettes, and the following:

Control material	REF
MRX Routine Normal Control	K5039

Solutions	REF
Deionised water for reconstitution e.g. MRX Laboratory Water	K5036

10 Quality control

To maintain consistent assay results, it is recommended that control plasmas are assayed at regular intervals. MRX Routine Normal Control (K5039) is recommended for MRX Thrombin Time. Each laboratory should establish a control range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice.

11 Results

The results are reported in seconds.

12 Expected values

Results from thrombin time measurements in plasma samples from 100 healthy blood donors using MRX Thrombin Time are presented below. The analysis was performed using an ACL Top 700 CTS instrument.

N	Mean TT	Max-min interval
100 samples	12.2 s	10.0 - 15.2 s

Due to inter-laboratory variation, each laboratory should set its own reference interval.

13 Limitations and interfering substances

The thrombin time may be prolonged in patients with an impaired polymerisation process, commonly observed in patients with hereditary dysfibrinogenemia.

The thrombin time may also be prolonged due to presence of fibrinogen degradation products and drugs, i.e., direct and indirect thrombin inhibitors.

The thrombin time is also affected by haemolysis and high lipid concentration.

MRX Thrombin Time is insensitive to the following substances on ACL Top 700 CTS instrument:

Interfering substance	Tolerance
Lipids	125 mg/dL
Bilirubin	40 mg/dL
Haemoglobin	250 mg/dL

Heparin sensitivity:

Determination of thrombin time in heparinised samples is dependent on the instrument and the kind of heparin used and may vary from laboratory to laboratory. Heparin tests were performed using a normal plasma pool spiked with un-fractionated heparin. Measurements were performed using an ACL Top 700 CTS instrument and the results show that MRX Thrombin Time is influenced by heparin:

Heparin (U/dL)	Lot A (s)	Lot B (s)
0	13.4	15.4
1	14.4	16.7
3	19.0	26.0/No clot
5 - 20	No clot	No clot

Dabigatran sensitivity:

Dabigatran tests were performed using a normal plasma pool spiked with dabigatran and analysed using an ACL Top 700 CTS instrument. The results show that MRX Thrombin Time is influenced by dabigatran:

Dabigatran (mg/dL)	Lot A (s)	Lot B (s)
0	13.7	15.1
0.001	32.2	38.4
0.002 – 0.04	No clot	No clot

Each laboratory is recommended to establish its own instrument-specific level where heparin and dabigatran is detected in citrated plasma.

14 Analytical performance characteristics

The following performance data was obtained with an ACL Top 700 CTS instrument. Performance will depend on the instrument used.

Precision:

Sample	Mean TT	Repeatability CV
Level 1	14.6 s	2.6%

15 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

16 Additional information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.
The instrument-specific application sheet is available from your local distributor.

17 References

- 1. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays. 6th ed. CLSI guideline H21. Clinical and Laboratory Standards Institute; 2024.

18 Definition of symbols

Manufacturer

Consult electronic instructions for use

nordicbiomarker.com/IFU

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Biological risks

Batch code

Contains biological material of animal origin

19 Revision history

Version	Changes to previous version
4.0	Section 4: Information on material of animal origin updated. Minor editorial revisions made.

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