HYPHEN BioMed

IVD
CE
2797

Instructions for Use [EN]

D-Dimer Control Low REF SC108K

D-Dimer Control High REF SC109K

For *In vitro* Diagnostic Use. *English, revision: 2026-10*

1. Intended use
For quality control of the assay D-Dimer Latex LRT. Intended to be used by professional laboratory personnel using analysers with turbidimetric detection.

2. Components
D-Dimer Control Low consists of:
10 x 1 mL lyophilised citrated human plasma enriched with D-Dimer.

D-Dimer Control High consists of:
10 x 1 mL lyophilised citrated human plasma enriched with D-Dimer.

3. Warnings and precautions
Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

D-Dimer Control Low and High contain material of human origin. Each donor has been tested by approved methods and found negative for the presence of HBsAg and anti-HIV I & II and anti-HCV. However, as no method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious.

D-Dimer Control Low contains material of animal origin and should be handled as potentially infectious.

4. Acceptable range
The D-Dimer concentration in each new lot of D-Dimer Control Low and High is determined with D-Dimer Latex LRT. Refer to the Certificate of Analysis for the lot-specific result and acceptable range. The acceptable range should be used as guidance only in the internal quality control procedure employed by the laboratory. Each laboratory is recommended to establish its own acceptable range.

5. Preparation
- Before opening, carefully tap the vial against a surface to collect the lyophilised material at the bottom.
- Add 1.00 mL deionised water. The water temperature should be 15 - 25 °C.
- Reseal the vial and let it stand for approximately 15 minutes at 15 - 25 °C.
- Gently mix by swirling or rotating until the content is completely reconstituted.

6. Storage and stability
Store at 2 - 8 °C. After reconstitution, stable in the closed original vial, provided no contamination occurs, for:

- 24 hours at 15 - 25 °C.
- 7 days at 2 - 8 °C.
- 4 weeks at -20 °C. Do not refreeze after thawing.

7. Material required but not provided
- Analyser with turbidimetric detection (refer to the Instructions for Use associated with D-Dimer Latex LRT for valid instruments).
- Deionised water for reconstitution.

Reagent	REF
D-Dimer Latex LRT	120011

8. Quality control
To maintain consistent assay results, it is recommended that controls in different levels are analysed together at regular intervals. Each laboratory is recommended to establish its own acceptable range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. A new acceptable range must be determined for each individual lot of controls. Recalibration is suggested at least when controls are not within the acceptable range and/or each time a new lot of reagent is used.

9. Summary of Safety and Performance
The Summary of Safety and Performance (SSP) is available in the European database on medical devices (Eudamed), where they are linked to these Basic UDI-DIs, 73500603203677 (SC108K), 73500603203779 (SC109K). The Eudamed public website is found on https://ec.europa.eu/tools/eudamed. In case Eudamed is not fully functional, the SSP is available from the manufacturer on request.

10. Reporting of incidents
Any serious incidents that occur in relation to these devices shall be reported to the manufacturer as well as the national competent authority in which the user is established.

11. Additional information
The e-IFU (other languages) and the safety data sheet are available at www.hyphen-biomed.com. A paper copy of these Instructions for Use is available on request. Contact your local distributor.

For customer support and instrument-specific application sheet, please contact your local distributor.

12. Definition of symbols

Manufacturer		www.hyphen-biomed.com	Consult electronic instructions for use
CE	CE mark		Use-by date
IVD	In vitro diagnostic medical device		Temperature limit
REF	Catalogue number		Biological risks
LOT	Batch code	BIO	Contains biological material of animal origin
UDI	Unique device identifier		Contains human blood or plasma derivatives

13. Revision history

Version	Changes to previous version
1.0	First version of document.

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Nordic Biomarker AB, Vildmannavägen 1, 903 47 Umeå, Suède