

POCT D-dimer Control

POCT D-dimer Control is a liquid plasma control material intended for use in External Quality Assurance (EQA) programs for D-dimer testing.

Developed for professional users and can be used on both automated analyzers and point-of-care (POC) systems.

The control is prepared from modified human plasma and can be supplied in various levels and volumes to meet specific laboratory requirements.



PRODUCT INFORMATION

Product form	Frozen liquid
Product origin	Human plasma
1 Month storage	2 - 8 °C
6 Months storage	-20 °C
Long term storage	≤ -70 °C

ORDERING INFORMATION

PRODUCT NUMBER	PRODUCT NAME	mg/L DDU	Volume
C3198	POCT D-dimer Control L1	<0.2	0.5 mL
C3199	POCT D-dimer Control L2	0.2 - 1.0	0.5 mL
C3200	POCT D-dimer Control L3	1.1 - 2.0	0.5 mL
C3201	POCT D-dimer Control L1 Bulk	<0.2	Custom vol.
C3202	POCT D-dimer Control L2 Bulk	0.2 - 1.0	Custom vol.
C3203	POCT D-dimer Control L3 Bulk	1.1 - 2.0	Custom vol.

Manufactured in Sweden

order@nordicbiomarker.com

www.nordicbiomarker.com