

POCT INR Control

POCT PT INR Control is a ready-to-use liquid plasma control intended for use in External Quality Assurance (EQA) programs for PT INR testing.

Developed for professional users and can be used on a wide range of POC systems. The control is supplied in five INR level ranges as single-dose 0.3 mL units in a convenient pipette bottle for immediate use after thawing.

Delivered as a frozen liquid formulation to ensure stable and consistent performance while minimizing preparation time and reducing the risk of handling errors.



PRODUCT INFORMATION

Product form	Frozen liquid
Product origin	Human plasma
1 Month storage	2 - 8 °C
6 Months storage	-20 °C
Long term storage	≤ -70 °C

ORDERING INFORMATION

PRODUCT NUMBER	PRODUCT NAME	INR	Volume
C3193	POCT INR Control L1	1.3 - 1.6	1 x 0.3 mL
C3194	POCT INR Control L2	1.7 - 2.0	1 x 0.3 mL
C3195	POCT INR Control L3	2.1 - 3.0	1 x 0.3 mL
C3196	POCT INR Control L4	3.1 - 4.0	1 x 0.3 mL
C3197	POCT INR Control L5	>4.0	1 x 0.3 mL

Manufactured in Sweden

order@nordicbiomarker.com

www.nordicbiomarker.com